

TABLE OF CONTENTS


1 Introduction	1
2 Installation	1
3 BUILD INDEX	2
4 FUSIONS DETECTION	
4.1 Prepare	3
4.2 Run arriba	3
4.3 OUTPUT FILE INTERPRETATION	4

1 Introduction

This pipeline could accurately and efficiently identify gene fusion from RNA-Seq data. It was developed for the use in a clinical research setting. Alignment is done by the ultrafast STAR aligner, and gene fusions detection tool arriba is based on STAR output result.

Workflow

Fusion detection with Arriba is based on the STAR aligner. It is an extension of the regular alignment workflow, which can be incorporated with few modifications. The addition of Arriba does not affect the normal alignments. The workflow yields fusion predictions as well as normal alignments that can be used for other downstream analyses such as expression quantification or variant calling. Like so, fusion detection incurs negligible computational overhead, since it adds only a few minutes of runtime to the regular alignment workflow.



2 Installation

Annotation Tools Recommendation – ANNOVAR



STAR

Spliced Transcripts Alignment to a Reference (STAR) is written in C and C++ and can be run as a standalone application on diverse hardware systems. It is freely available to personal, academic and non-profit use only. You cannot redistribute ANNOVAR to other users including lab members. Download the latest package from the STAR <u>releases page</u> as shown:

Example:

```
wget https://github.com/alexdobin/STAR/archive/refs/tags/2.7.10a.tar.gz
tar -xzf 2.7.10a.tar.gz
cd STAR-2.7.10a/source
make STAR # Compile under Linux
```

After compilation, add STAR to environment variable. #remember to change the path to the absolute path to STAR-2.7.10a/source in your server.

Example:

```
export PATH=/mnt/sata/Janet_Analysis/GeneFusion/STAR-2.7.10a/source/:$PATH
```

ARRIBA

Arriba has only a single prerequisite: STAR (version >=2.7.10a recommended). Download the latest tarball from the arriba <u>releases page</u> as shown.

Example:

```
wget https://github.com/suhrig/arriba/releases/download/v2.2.1/arriba_v2.2.1.tar.gz
tar -xzf arriba_v2.2.1.tar.gz
cd arriba_v2.2.1 && make
```

Arriba requires an assembly in FastA format, gene annotation in GTF format, and a STAR index built from the two. So STAR index need to be built before running arriba.

3 BUILD INDEX

If you do not already have the files and a STAR index, you can use the script <code>downLoad_references.sh</code>. It downloads the files to the current working directory and builds a STAR index. GENCODE annotation is recommended over RefSeq due to more comprehensive annotation of immunoglobulin/T-cell receptor loci and splice sites, which improves sensitivity.

Run the script without arguments to see a list of available files.

Annotation Tools Recommendation – ANNOVAR



```
(base) [ztron@MegaBOLT Workstation arriba v2.2.1]$ ./download references.sh hg19+RefSeq hg19
Usage: download_references.sh ASSEMBLY+ANNOTATION
Available assemblies and annotations:
GRCh37+ENSEMBL87
GRCh37+GENCODE19
GRCh37+RefSeq
GRCh37viral+ENSEMBL87
GRCh37viral+GENCODE19
GRCh37viral+RefSeg
GRCh38+ENSEMBL93
GRCh38+GENCODE28
GRCh38+RefSeg
GRCh38viral+ENSEMBL93
GRCh38viral+GENCODE28
GRCh38viral+RefSeg
GRCm38+GENCODEM25
GRCm38+RefSeq
GRCm38viral+GENCODEM25
GRCm38viral+RefSeg
GRCm39+GENCODEM26
GRCm39+RefSeq
GRCm39viral+GENCODEM26
GRCm39viral+RefSeq
hg19+ENSEMBL87
hg19+GENCODE19
hg19+RefSeq
hg19viral+ENSEMBL87
hg19viral+GENCODE19
hg19viral+RefSeq
hg38+ENSEMBL93
```

Download reference files and build index:

```
./download references.sh $ASSEMBLIES+ANNOTATIONS
```

Example:

```
./download_references.sh hg19+GENCODE19
```

Output:

4 FUSIONS DETECTION

4.1 Prepare

Run the script without arguments to check usage.

Annotation Tools Recommendation - ANNOVAR



(base) [ztron@MegaBOLT_Workstation arriba_v2.2.1]\$./run_arriba.sh
Usage: run_arriba.sh STAR_genomeDir/ annotation.gtf assembly.fa blacklist.tsv known_fusions.tsv protein_domains.gff3 threads read1.fastq.gz [read2.fastq.gz]

blacklist.tsv, known_fusions.tsv and protein_domains.gff3 can be found under directory arriba_v2.2.1/database/. Choose the corresponding files according to the coordinates (eg. hg19/hs37d5/GRCh37 shares the same chromosome coordinates, hg38/GRCh38 shares the the same chromosome coordinates and vice versa). Unzip files needed using command line:

```
gunzip filename.gz
```

Example:

```
gunzip database/blacklist_hg19_hs37d5_GRCh37_v2.2.1.tsv.gz
gunzip database/known_fusions_hg19_hs37d5_GRCh37_v2.2.1.tsv.gz
```

Output:

```
(base) [ztron@MegaBOLT_Workstation arriba_v2.2.1]$ ls database/
blacklist_hg19 hs37d5 [GRCh37_v2.2.1.tsv optobands_hg38 [GRCh38_v2.2.1.tsv optobands_hg38] [GRCh38_v2.2.1.tsv optobands_hg38] [GRCh38_v2.2.1.tsv optobands_hg38] [GRCh38_v2.2.1.tsv optobands_hg38] [GRCh38_v2.2.1.tsv optobands_hg19_hs37d5 [GRCh37_v2.2.1.tsv optobands_hg19_hs37d5 [GRCh37_v2.2.1.tsv optobands_hg19_hs37d5 [GRCh37_v2.2.1.tsv optobands_hg19_hs37d5 [GRCh38_v2.2.1.tsv] optobands_hg19_hs37d5 [GRCh38_v2.2.1.tsv] optobands_hg19_hs37d5 [GRCh37_v2.2.1.tsv] optoba
```

If you use other assemblies whose coordinates are incompatible with hg19/hs37d5/GRCh37 or hg38/GRCh38 or mm10/GRCm38 or mm39/GRCm39, then the coordinates in the blacklist will not match and the predictions will contain many false positives.

4.2 Run arriba

Run arriba by demo script run arriba.sh

Example:

(Run the demo script with 8 threads)

```
sh run_arriba.sh STAR_index_hg19_RefSeq_hg19/ RefSeq_hg19.gtf hg19.fa \
database/blacklist_hg19_hs37d5_GRCh37_v2.2.1.tsv \
database/known_fusions_hg19_hs37d5_GRCh37_v2.2.1.tsv \
database/protein_domains_hg19_hs37d5_GRCh37_v2.2.1.gff3 8 read1.fq.gz read2.fq.gz
```

4.3 Output File Interpretation

See: https://arriba.readthedocs.io/en/latest/output-files/